Keep Your Instruments Up and Running with CarboPrep **Plus SPE Cartridges**

- Maximum Contaminant Removal
- Less Instrument Downtime

Resprep

Higher Sample Throughput

Pure Chromatography

RESTÊK

Resprep

CarboPrep Plus

25845

Resprep

CarboPrep Plus

CarboPrep Plus

25845

www.restek.com


Four Words No Environmental Lab Manager Wants to Hear: "The Instrument is Down."

If you are in charge of an environmental lab, you know that instruments that are down for unscheduled maintenance due to inlet or column contamination are not only costing money instead of making it, but are also jeopardizing your ability to meet your customers' deadlines and performance requirements.

Chances are, you are also in charge of the lab that prepares samples for analysis, so you appreciate the value in making that process simple, reproducible, and effective, mitigating the chance for error and providing a clean extract that won't shut an instrument down in the analytical lab because of inlet-fouling matrix components that made their way through the sample cleanup process.

Unfortunately, for the analysis of organochlorine pesticides in soil and water, the industry standard sample preparation product—activated Florisil SPE cartridges—just isn't up to the job of removing the worst of those sample matrix contaminants that condense inside inlet liners or coat the inlet side of your analytical column. Florisil is great for stopping polar compounds, but those usually aren't the compounds that cause instruments to fail their calibration checks, requiring unscheduled maintenance, delaying processing new samples, and adding the expense of replacing consumables ahead of schedule.

If your samples may contain those relatively high molecular weight sample matrix compounds that are the bane of injection ports and sensitive analytes, you'll need something different. Fortunately, new Resprep CarboPrep Plus SPE cartridges outperform Florisil cartridges in removing these troublesome compounds, protecting your system and keeping samples running through your lab efficiently.






An Easy Change for the Sample Prep Lab

To remove the matrix components that contaminate your inlets and degrade your analytes of interest, you need a sample prep product with a different selectivity than Florisil that still allows all the target analytes to elute. But, chances are, your lab is set up to run a standard Florisil cleanup and switching to a new process can be a logistical challenge.

To make the switch to another SPE product easier, Restek's CarboPrep Plus cartridges use the same equipment, solvents, and solvent volumes that work for Florisil, but the specially treated carbon in the CarboPrep Plus tubes will do a much better job of removing those inlet-contaminating matrix components that may be present in your customers' samples.

While the CarboPrep Plus SPE cartridges have different selectivity than activated Florisil cartridges, they will still elute the compounds of interest and remove contaminating compounds using the same procedure. Change doesn't get easier than that.





Rock-Solid Reliability Also Benefits the Analytical Lab


A proprietary treatment renders the carbon in CarboPrep Plus SPE cartridges clean and consistent, allowing you to standardize on an extraction protocol and not have to optimize elution volumes for every new batch of cartridges. This simplifies workflow in the sample prep lab and provides confidence that the next cartridge will perform exactly the same as the last. The analytical lab also benefits from a CarboPrep Plus cleanup because less downtime is required for maintenance, and high recoveries and low background levels are consistently obtained (Table I). Less maintenance and excellent analytical performance mean higher sample throughput for your lab.

> % RSD Compound Average % Recovery* 1.3 TCMX 99.8 hexachlorobenzene 97.6 2.3 alpha-BHC 99.2 1.3 1.4 gamma-BHC 98.1 beta-BHC 98.7 1.4 delta-BHC 1.7 98.7 heptachlor 101.2 1.7 aldrin 97.1 1.7 heptachlor epoxide 98.6 1.6 gamma-chlordane 98.2 1.8 alpha-chlordane 96.9 2.2 4,4'-DDE 96.8 1.6 endosulfan I 98.0 1.4 dieldrin 97.8 1.4 endrin 98.4 1.6 4,4'-DDD 1.4 97.4 endosulfan II 94.6 5.5 4,4'-DDT 98.2 1.2 endrin aldehyde 96.1 1.2 100.2 2.1 methoxychlor 2.2 endosulfan sulfate 97.1 endrin ketone 98.1 0.8 DCB 97.8 2.2 * n = the analysis of 6 cartridges each from 3 different lots for a total of 18 data points


Table I: Accurate, consistent recoveries for organochlorine pesticides and surrogate compounds (TCMX and DCB) make switching to CarboPrep Plus cartridges an easy way to improve analytical performance.



Figure 3: Lot-to-lot consistency ensures that all the compounds of interest will elute within the 10 mL collection volume used for typical cleanup procedures, ensuring excellent recoveries time and again. In this figure, you can see a subset of organochlorine pesticides that are retained to different degrees on CarboPrep Plus material, and even hexachlorobenzene, the most highly retained compound, showed consistently high recovery at 97.6% (2.3% RSD, Table I). Note that all compounds listed in Table I elute from CarboPrep Plus cartridges prior to tetrachlorometaxylene (TCMX), so elution volumes can be reduced to save time for labs not analyzing hexachlorobenzene.







Make the Simple Switch to CarboPrep Plus SPE Cartridges and Improve Lab Productivity

If you have trouble with samples contaminating your GC inlets and columns, resulting in failing calibration checks for your organochlorine pesticides, unscheduled instrument maintenance, rework, delays, and unexpected expenses, then consider adding CarboPrep Plus cartridges to your toolbox of sample preparation options. While Florisil is excellent at removing polar matrix contaminants, the matrix contaminating compounds that foul inlets and columns pass through and end up in your instruments. To get consistent recoveries, while mitigating the risk that a dirty sample will take an instrument down, choose CarboPrep Plus SPE cartridges for your organochlorine pesticides workflow.



Resprep CarboPrep Plus SPE Cartridges

- Designed specifically for the cleanup of sample extracts for organochlorine pesticides analysis.
- Excellent alternative to Florisil products, especially for the removal of nonvolatile matrix components that contaminate GC inlets and columns.
- Proprietary treatment renders the carbon consistent and clean, ensuring the same selectivity tube to tube, and lot to lot with no interfering background.
- Uses the same hardware, solvents, and solvent volumes as traditional Florisil cleanup, so switching is simple.

Description	Volume	qty.	cat.#
Resprep CarboPrep Plus SPE Cartridges	3 mL, 95 mg	30-pk.	25845



Related Products

Resprep 12- or 24-Port SPE Manifolds

Complete Resprep SPE Manifold, 12-Port (Includes: glass basin with built-in vacuum regulator (1); polypropylene lid with 12 individual control valves (1); 12-position collection rack (1); PTFE sample guides (12); and waste container (1).)kit26077Complete Resprep SPE Manifold, 24-Port (Includes: glass basin with built-in vacuum regulator (1); polypropylene lid with 24 individual control valves (1); 24-position collection rack (1); and PTFE sample guides (24).)kit26080Resprep Manifold Replacement Parts </th <th>Description</th> <th>qty.</th> <th>cat.#</th>	Description	qty.	cat.#
(Includes: glass basin with built-in vacuum regulator (1); polypropylene lid with 24 individual control valves (1); 24-position collection rack (1); and PTFE sample guides (24).)kit26080Resprep Manifold Replacement Parts dy. cat.#Replacement Waste Container, 12-Portea.24014Replacement Vacuum Valve and Gauge Assemblyea.24008Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Portea.25991 Collection Rack Collection Rack, 12-Portea.26079Collection Rack, 12-Portea.26082Manifold LidReplacement Manifold Lid (sample guides not included), 12-Portea.26082Manifold LidReplacement PartsValves, 12 or 24-Port2-pk.26083Valves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.26084	(Includes: glass basin with built-in vacuum regulator (1); polypropylene lid with 12 individual control valves (1); 12-position collection rack (1); PTFE sample guides (12); and waste	kit	26077
Descriptionqty.cat.#Replacement Waste Container, 12-Portea.24014Replacement Vacuum Valve and Gauge Assemblyea.24008Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Portea.25991Collection RackCollection Rack, 12-Portea.26079Collection Rack, 24-Portea.26082Manifold Lidmanifold Lid (sample guides not included), 12-Portea.26082Manifold Lid Replacement Manifold Lid (sample guides not included), 24-Portea.26081Valves, 12 or 24-Port2-pk.2608326083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.2608426084	(Includes: glass basin with built-in vacuum regulator (1); polypropylene lid with 24 individual	kit	26080
Replacement Waste Container, 12-Portea.24014Replacement Vacuum Valve and Gauge Assemblyea.24008Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Portea.25991Collection RackCollection Rack, 12-Portea.26079Collection Rack, 24-Portea.26082Manifold LidReplacement Manifold Lid (sample guides not included), 12-Portea.26073Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement PartsValves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.26084	Resprep Manifold Replacement Parts		
Replacement Vacuum Valve and Gauge Assemblyea.24008Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Portea.25991Collection RackECollection Rack, 12-Portea.26079Collection Rack, 24-Portea.26082Manifold LidReplacement Manifold Lid (sample guides not included), 12-Portea.26078Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement Parts2-pk.26083Valves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.26084	Description	qty.	cat.#
Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Portea.25991Collection RackECollection Rack, 12-Portea.26079Collection Rack, 24-Portea.26082Manifold LidReplacement Manifold Lid (sample guides not included), 12-Portea.26078Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement Parts2-pk.26083Valves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.26084	Replacement Waste Container, 12-Port	ea.	24014
Collection RackCollection Rack, 12-Portea.26079Collection Rack, 24-Portea.26082Manifold LidReplacement Manifold Lid (sample guides not included), 12-Portea.26078Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement PartsValves, 12 or 24-Port2-pk.26083Valves, 12 or 24-Port2-pk.2608326130Needle Attachment48-pk.2613126084	Replacement Vacuum Valve and Gauge Assembly	ea.	24008
Collection Rack, 12-Portea.26079Collection Rack, 24-Portea.26082Manifold LidEan Collection Rack, 24-Portea.26078Replacement Manifold Lid (sample guides not included), 12-Portea.26078Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement PartsValves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.2608426084	Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Port	ea.	25991
Collection Rack, 24-Portea.26082Manifold LidExplacement Manifold Lid (sample guides not included), 12-Portea.26078Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement PartsValves, 12 or 24-Port2-pk.26083Valves, 12 or 24-Port2-pk.26030Needle Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.26084	Collection Rack		
Manifold LidReplacement Manifold Lid (sample guides not included), 12-Portea.26078Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement PartsValves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.2608426084	Collection Rack, 12-Port	ea.	26079
Replacement Manifold Lid (sample guides not included), 12-Portea.26078Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement PartsValves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.26084	Collection Rack, 24-Port	ea.	26082
Replacement Manifold Lid (sample guides not included), 24-Portea.26081Manifold Lid Replacement PartsValves, 12 or 24-Port2-pk.26083Valve Attachment48-pk.26130Needle Attachment48-pk.26131Sample Guide Needles, 12 or 24-Port12-pk.26084	Manifold Lid		
Manifold Lid Replacement Parts Valves, 12 or 24-Port 2-pk. 26083 Valve Attachment 48-pk. 26130 Needle Attachment 48-pk. 26131 Sample Guide Needles, 12 or 24-Port 12-pk. 26084	Replacement Manifold Lid (sample guides not included), 12-Port	ea.	26078
Valves, 12 or 24-Port 2-pk. 26083 Valve Attachment 48-pk. 26130 Needle Attachment 48-pk. 26131 Sample Guide Needles, 12 or 24-Port 12-pk. 26084	Replacement Manifold Lid (sample guides not included), 24-Port	ea.	26081
Valve Attachment 48-pk. 26130 Needle Attachment 48-pk. 26131 Sample Guide Needles, 12 or 24-Port 12-pk. 26084	Manifold Lid Replacement Parts		
Needle Attachment 48-pk. 26131 Sample Guide Needles, 12 or 24-Port 12-pk. 26084	Valves, 12 or 24-Port	2-pk.	26083
Sample Guide Needles, 12 or 24-Port 12-pk. 26084	Valve Attachment	48-pk.	26130
	Needle Attachment	48-pk.	26131
Gasket, 12-port 2-pk. 24011	Sample Guide Needles, 12 or 24-Port	12-pk.	26084
	Gasket, 12-port	2-pk.	24011



Rtx-CLPesticides Columns (fused silica)

Description	temp. limits	qty.	cat.#
Rtx-CLPesticides Column 30 m, 0.25 mm ID, 0.25 µm	-60 to 320/340 °C	ea.	11123
Rtx-CLPesticides Column 30 m, 0.32 mm ID, 0.32 µm	-60 to 320/340 °C	ea.	11141
Rtx-CLPesticides Column 30 m, 0.53 mm ID, 0.50 µm	-60 to 300/320 °C	ea.	11140

Rtx-CLPesticides2 Columns (fused silica)

Description	temp. limits	qty.	cat.#
Rtx-CLPesticides2 Column 30 m, 0.25 mm ID, 0.20 µm	-60 to 320/340 °C	ea.	11323
Rtx-CLPesticides2 Column 30 m, 0.32 mm ID, 0.25 μm	-60 to 320/340 °C	ea.	11324
Rtx-CLPesticides2 Column 30 m, 0.53 mm ID, 0.42 μm	-60 to 300/320 °C	ea.	11340

2.0 mL, 9 mm Short-Cap, Screw-Thread Vials (vial only)

Fit all 2.0 mL, 12 x 32 mm, crimp-top vial-based autosamplers.


		100-pk.	1,000-pk.
Description	Color	cat.#	cat.#
Short-Cap Vial w/White Graduated Marking Spot	Clear	21140	21141
Short-Cap Vial w/White Graduated Marking Spot	Amber	21142	21143

2.0 mL, 11 mm Crimp-Top Vials, 12 x 32 mm

Choose with or without graduated marking spots.

Description	Color	Material	100-pk. cat.#	1,000-pk. cat.#
Vial w/White Graduated Marking Spot	Clear	Glass	24383	24384
Vial w/White Graduated Marking Spot	Amber	Glass	24385	24386

Other vials available! Find the right one at www.restek.com/vials



Also available in money-saving kits!





2.0 mL, 9 mm Short-Cap, Screw-Vial Closures (Polypropylene, preassembled)

			100-pk.	1,000-pk.
Cap Color	Septa Material	Туре	cat.#	cat.#
Blue	PTFE/Silicone	Screw-Thread	24485	24486

2.0 mL, 11 mm Aluminum Crimp Seals With Septa

		100-pk.	1,000-pk.
Seal Color	Septa Material	cat.#	cat.#
Blue	PTFE/Silicone	24361	24362

2.0 mL, 9 mm Screw-Thread Vial Convenience Kits

(Vials, Caps, & Septa)

Vials packaged in a clear-lid tray. Caps with septa packaged in a plastic bag.

	100-pk.	1,000-pk.
Description	cat.#	cat.#
Clear Vial w/Graduated Marking Spot, 9 mm, Blue Cap PTFE/Silicone Liner	26592	26593
Amber Vial w/Graduated Marking Spot, 9 mm, Blue Cap PTFE/Silicone Liner	26596	26597

2.0 mL, 11 mm Crimp Vial Convenience Kits

(Vials, Caps, & Septa)

Vials packaged in a clear-lid tray. Caps with septa packaged in a plastic bag.

	100-pk.	1,000-pk.
Description	cat.#	cat.#
Clear 2.0 mL Vial, Untreated, Silver Seal, PTFE/Silicone Septa	24646	24647
Amber 2.0 mL Vial, Untreated, Silver Seal, PTFE/Silicone Septa	24648	24649

Topaz 4.0 mm ID Single Taper Inlet Liner w/ Wool

for Agilent GCs equipped with split/splitless inlets

RESTEK		
ID x OD x Length	qty.	cat.#
Single Taper, Premium Deactivation, Borosilicat	e Glass with Quartz Wool	
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23303

Topaz 3.5 mm ID Single Taper Inlet Liner w/ Wool

for Shimadzu 17A, 2010, 2014, and 2030 GCs equipped with split/ splitless inlets

RESTEK		
ID x OD x Length	qty.	cat.#
Single Taper, Premium Deactivation, Borosilicate Gl	ass with Quartz Wool	
3.5 mm x 5.0 mm x 95 mm	5-pk.	23336

Topaz 4.0 mm ID Single Taper Inlet Liner w/ Wool

for Thermo TRACE 1300/1310 GCs equipped with SSL inlets

RESTEK		
ID x OD x Length	qty.	cat.#
Single Taper, Premium Deactivation, Boros	ilicate Glass with Quartz V	Nool
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23447

Organochlorine Pesticide Mix AB #2

- j	
(20 components) Aldrin (309-00-2), 8 μg/mL α-BHC (319-84-6), 8 μg/mL β-BHC (319-85-7), 8 μg/mL β-BHC (319-86-8), 8 μg/mL γ-BHC (Lindane) (58-89-9), 8 μg/mL <i>cis</i> -Chlordane (5103-71-9), 8 μg/mL <i>trans</i> -Chlordane (5103-74-2), 8 μg/mL 4,4'-DDD (72-54-8), 16 μg/mL 4,4'-DDT (50-29-3), 16 μg/mL 4,4'-DDT (50-29-3), 16 μg/mL Endosulfan I (959-98-8), 8 μg/mL Endosulfan II (33213-65-9), 16 μg/mL Endosulfan sulfate (1031-07-8), 16 μg/mL Endosulfan sulfate (1031-07-8), 16 μg/mL	
Endosulfan sulfate (1031-07-8), 16 µg/mL	
Endrin (72-20-8), 16 µg/mL Endrin aldehyde (7421-93-4), 16 µg/mL	
Endrin ketone (53494-70-5), 16 µg/mL	
Heptachlor (76-44-8), 8 µg/mL Heptachlor epoxide (isomer B) (1024-57-3), 8 µg/mL Methoxychlor (72-43-5), 80 µg/mL	
In hexane:toluene (1:1), 1 mL/ampul	cat.# 32292 (ea.)

2,4,5,6-Tetrachloro-m-xylene

2,4,5,6-Tetrachloro-*m*-xylene (877-09-8)

200 µg/mL in acetone, 1 mL/ampul	cat.# 32027 (ea.)
200 µg/mL in acetone, 5 mL/ampul	cat.# 32028 (ea.)

Decachlorobiphenyl, 508A

Decachlorobiphenyl (BZ #209) (2051-24-3)

200 µg/mL in acetone, 1 mL/ampul	cat.# 32029 (ea.)
200 µg/mL in acetone, 5 mL/ampul	cat.# 32030 (ea.)
10 µg/mL in isooctane, 1 mL/ampul	cat.# 32289 (ea.)

Hexachlorobenzene

Hexachlorobenzene (118-74-1)

1,000 µg/mL in acetone, 1 mL/ampul	cat.# 32231 (ea.)

Certified reference materials (CRMs) manufactured and QC-tested in ISO-accredited labs satisfy your ISO requirements.

Use our Rtx-CLPesticides/Rtx-CLPesticides2 capillary column pair to analyze all twenty organochlorine pesticides and their surrogates simultaneously.



Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

 ${\scriptstyle \circledcirc}$ 2018 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# EVSS2976-UNV